

# Fluorescence Spectroscopy

03/12/02: Marcia Levitus

- absorption of light.
- radiative and non-radiative processes in molecules.
- properties of fluorescence.
- experimental determinations.
- fluorescence of macromolecules
  - > intrinsic and extrinsic fluorophores
  - > solvent effects
  - > quenching

03/14/02: Carlos Bustamante

Energy transfer and fluorescence polarization

03/19/02: Marcia Levitus

Single molecule fluorescence



# Interaction of light with matter

The interaction of radiation with matter can cause transitions between the energy levels of the atoms or molecules. The possible energy of an electronic orbital, molecular vibration, or rotation is restricted to a well-defined energy.


Type of Radiation	Wavelength Range	Type of Transition
gamma-rays	<1 pm	nuclear
X-rays	1 nm-1 pm	inner electron
ultraviolet	400 nm-1 nm	outer electron
visible	750 nm-400 nm	outer electron
near-infrared	2.5 $\mu$ m-750 nm	outer electron molecular vibrations
infrared	25 $\mu$ m-2.5 $\mu$ m	molecular vibrations
microwaves	1 mm-25 $\mu$ m	molecular rotations, electron spin flips
radio waves	>1 mm	nuclear spin flips

$$k_B T (300K) = 25 \text{ meV} \rightarrow \approx 50 \mu\text{m}$$

————— Electronic levels  
 ————— Vibrational levels



# Absorption of light



$$A = -\log(T) = -\log\left(\frac{I}{I_0}\right)$$

$$I_0 = I + I_{abs} \quad (\text{no reflected light})$$

$$I_{abs} = I_0(1 - 10^{-A})$$



$$A = \epsilon \cdot l \cdot C$$

$\epsilon$ : extinction coefficient ( $\text{M}^{-1} \cdot \text{cm}^{-1}$ )

$l$ : path length (cm)




$C$ : concentration (M)

A=Absorbance, T= Transmittance




# Jablonski diagram


Internal conversion: in condensed phases molecules rapidly relax to the ground vibrational state of the first excited electronic state




## Examples of fluorophores






**YOYO-1 bound to DNA.**



**fluorescein**



**ethidium bromide  
bound to DNA.**



## Properties of Fluorescence Emission:

1- The absorption spectrum reflects the vibrational levels of the excited electronic state. The fluorescence spectrum reflects the vibrational levels of the ground state.

2- Stokes shift:

$$\lambda_{\text{emission}} > \lambda_{\text{absorption}}$$


(even for the  $0 \rightarrow 0$  transition)

3- The spectrum of the emitted light is invariant with the excitation wavelength. because of the fast relaxation to the lowest vibrational level of  $S_1$


Exception: azulene, it can emit from both  $S_1$  and  $S_2$ .

4- The mirror image rule.


The fluorescence emission often appears as the mirror image of the absorption spectrum.



The Franck-Condon principle: transitions are vertical in both absorption and emission




Nuclear Displacement



Nuclear Displacement


The Franck-Condon factor is the same for absorption and fluorescence



Exceptions:


- very long lived  $S_1$  state : emission occurs from a different geometry.
- Reactions from the excited state.

## Fluorescence lifetime and quantum yield



In general the life-time of the excited state is shorter than  $\tau_r$  due to other deactivation pathways from  $S_1$ .





$k_r$  = fluorescence radiative rate

$k_{ic}$  = internal conversion (collision with solvent + internal vibrational modes)  
 $\rightarrow$  increases with T

$k_Q$  = quenching (collision with solute molecules, such as  $O_2$ )

$k_T$  = energy transfer (FRET)

$k_x$  = reaction from excited state

$k_{ISC}$  = intersystem crossing (forbidden)

$$k = (k_r + k_{IC} + k_{ISC} + k_x + k_Q[Q] + k_{FRET} \dots)$$

$$-\frac{d[S_1]}{dt} = k[S_1] \rightarrow [S_1](t) = [S_1](0)e^{-t/\tau_f}$$

$$\tau_f = (k_r + k_{IC} + k_{ISC} + k_x + k_Q[Q] + k_{FRET} \dots)^{-1}$$


$$\phi = \frac{\text{photons emitted}}{\text{photons absorbed}} = \frac{k_r}{k} = \tau_f \cdot k_r = \frac{\tau_f}{\tau_r}$$

$\phi$  = fluorescence quantum yield;


$\tau_f$  = fluorescence life-time

$\tau_r$  = radiative fluorescence life-time


## Effect of structure and environmental rigidity



$\phi$  (25°C, organic solvent) = 0.05  
 $\phi$  (25°C, glycerol) = 0.15  
 $\phi$  (77K, organic solvent) = 0.75




$\phi(25^\circ\text{C}) = 1$   
 $\phi(77\text{K}) = 1$



— □  
 — □ 77K  
 — □ 25°C

$$\phi = \frac{\text{photons emitted}}{\text{photons absorbed}} = \frac{k_r}{k} = \tau_f \cdot k_r = \frac{\tau_f}{\tau_r}$$


## Measurement of fluorescence spectra




**M1: excitation monochromator**

**M2: emission monochromator**


**Absorption spectrum**



**Fl. Emission spectrum**



**Fl. Excitation spectrum**



## Measurement of fluorescence quantum yields

fraction of intensity emitted at that particular wavelength


fraction of total fluorescence that is detected

$$I_f(\lambda_{em}) = I_{Abs}(\lambda_{ex}) \cdot \phi_f \cdot f(\lambda_{em}) \cdot K$$

fluorescence quantum yield

absorbed intensity at  $\lambda_{ex}$

measured intensity of fluorescence at  $\lambda_{em}$



$$I_{abs}(\lambda_{ex}) = I_0(\lambda_{ex}) (1 - 10^{-A(\lambda_{ex})})$$

If  $A \rightarrow 0$   $I_{abs}(\lambda_{ex}) \cong 2.303 \cdot A(\lambda_{ex}) \cdot I_0(\lambda_{ex})$

$$\square_{\lambda_{em}} I_f = I_{abs}(\lambda_{ex}) \cdot \phi_f \cdot K \cong 2.303 \cdot A(\lambda_{ex}) \cdot I_0(\lambda_{ex}) \cdot \phi_f \cdot K$$

**If we measure the sample and a standard under the same experimental conditions, keeping  $\lambda_{ex}$  constant:**

Standards:

Quinine sulfate in  $H_2SO_4$  1N:  $\phi_f = 0.55$

Fluorescein in NaOH 0.1N:  $\phi_f = 0.93$


$$\frac{\square_{\lambda_{em}} I_f^{\text{sample}}}{\square_{\lambda_{em}} I_f^{\text{standard}}} = \frac{A^{\text{sample}}(\lambda_{ex})}{A^{\text{standard}}(\lambda_{ex})} \cdot \frac{\phi_f^{\text{sample}}}{\phi_f^{\text{standard}}}$$

Important : the index of refraction of the two solvents (sample and standard) must be the same



## Measurement of fluorescence lifetimes

### Time correlated single photon counting:

The TCSPC measurement relies on the concept that the probability distribution for emission of a single photon after an excitation yields the actual intensity against time distribution of all the photons emitted as a result of the excitation. By sampling the single photon emission after a large number of excitation flashes, the experiment constructs this probability distribution.



pulsed source



## Intrinsic fluorescence of macromolecules

Substance	Conditions	Absorption		Fluorescence <sup>1</sup>		Sensitivity $\rightarrow \epsilon_{\max} \cdot \phi_f$	
		$\lambda_{\max}$ (nm)	$\epsilon_{\max}$ $\times 10^{-3}$	$\lambda_{\max}$ (nm)	$\phi_f$	$\tau_f$ (nsec)	$\epsilon_{\max} \phi_f$ $\times 10^{-3}$
→ Tryptophan	H <sub>2</sub> O, pH 7	280	5.6	348	0.20	2.6	11.
Tyrosine	H <sub>2</sub> O, pH 7	274	1.4	303	0.14	3.6	2.0
Phenylalanine	H <sub>2</sub> O, pH 7	257	0.2	282	0.04	6.4	0.08
Y base	Yeast tRNA <sup>Phe</sup>	320	1.3	460	0.07	6.3	0.91
Adenine	H <sub>2</sub> O, pH 7	260	13.4	321	$2.6 \times 10^{-4}$	<0.02	0.032
Guanine	H <sub>2</sub> O, pH 7	275	8.1	329	$3.0 \times 10^{-4}$	<0.02	0.024
Cytosine	H <sub>2</sub> O, pH 7	267	6.1	313	$0.8 \times 10^{-4}$	<0.02	0.005
Uracil	H <sub>2</sub> O, pH 7	260	9.5	308	$0.4 \times 10^{-4}$	<0.02	0.004
NADH	H <sub>2</sub> O, pH 7	340	6.2	470	0.019	0.40	1.2

Emission maxima and lifetimes of single **tryptophan** proteins



Usually determines the fluorescence of proteins.  
Tyr often is quenched by nearby tryptophanes


Protein	$\lambda_{\max}$ (nm)	$\tau_0$ (nsec)
Azurin	308	4.0
RNaseT1	324	3.5
HAS	342	6.0
Nuclease	334	5.0
Monellin	342	2.6
Glucagon	352	2.8
ACTH	352	3.1

## Protein fluorescence

Tryptophan, the dominant intrinsic fluorophore, is generally present at about 1 mol% in proteins. A protein may possess just one or a few Trp residues, which facilitates interpretation of the spectral data.


Trp is very sensitive to its local environment. It is possible to see changes in emission spectra in response to conformational changes, subunit association, substrate binding, denaturation, and anything that affects the local environment surrounding the indole ring. Also, Trp appears to be uniquely sensitive to collisional quenching, either by externally added quenchers, or by nearby groups in the protein.

Trp fluorescence can be selectively excited at 295-305 nm. (to avoid excitation of Tyr)



## Protein fluorescence: solvent effects.

General solvent effects:



Lippert equation:

$$(\bar{\nu}_A - \bar{\nu}_F) = \frac{2}{hC} \left( \frac{\epsilon - 1}{2\epsilon + 1} - \frac{n^2 - 1}{2n^2 + 1} \right) \frac{(\mu_E - \mu_G)^2}{a^3} + \text{constant}$$

Trp displays complex spectral properties due to the presence of two nearly isoenergetic excited states. These electronic transitions display distinct absorption, and emission spectra, and are differently sensitive to solvent polarity





Figure 16.5. Emission spectra of indole in cyclohexane, ethanol, and their mixtures at 20 °C. From Ref. 23.





Emission spectra of *Pseudomonas fluorescens* azurin Pfl.  
For 275-nm excitation, a peak is observed due to the tyrosine residue(s)

The position and structure of the fluorescence suggests that the indole residue is located in a completely nonpolar region of the protein. These results agree with X-ray studies, which show that the indole group is located in the hydrophobic core of the protein.




In the presence of a denaturing agent, the TrpP emission loses its structure and shifts to 351 nm, characteristic of a fully exposed Trp residue.

Changes in emission spectra can be used to follow protein unfolding

## Time-resolved protein fluorescence:

Resolution of the contributions of individual tryptophan residues in multi-tryptophan proteins.

$$I(\lambda, t) = \sum_i \alpha_i(\lambda) \exp(-t/\tau_i)$$



Collisional quenching of proteins has been extensively utilized to determine the extent of Trp exposure to the aqueous phase.

Collisional quenching requires the fluorophore and quencher to be in molecular contact. A water-soluble quencher will not easily penetrate the protein matrix.

$k_Q$  reflects the efficiency of quenching, or the accessibility of the fluorophores to the quencher.

Diffusion-controlled quenching typically results in values of  $k_Q$  near  $10^{10} \text{ M}^{-1}\text{s}^{-1}$ .








Figure 16.26. Collisional quenching of surface-accessible ( $W_1$ ) and buried ( $W_2$ ) tryptophan residues in proteins.

$$\frac{\phi}{\phi_0} = \frac{k_r / (k_r + k_{IC} + k_Q[Q])}{k_r / (k_r + k_{IC})}$$

without quencher

$$\frac{\phi_0}{\phi} = 1 + \frac{k_Q}{k_r + k_{IC}}[Q] = 1 + \tau_0 k_Q[Q] = \frac{I_0}{I} = \frac{\tau_0}{\tau}$$



Apoazurin Ade has two Trp residues, one surface residue and one buried residue.


In the presence of iodide (collisional quencher), the spectrum resembles the structured emission seen in non-polar environments.

The spectrum of the quenched Trp can be seen from the difference spectrum (dashed line), and is characteristic of an exposed residue in a partially hydrophobic environment.

**Green fluorescent protein** (abbreviated GFP) is isolated from the Pacific jellyfish *Aequorea victoria* and now plays central roles in biochemistry and cell biology due to its widespread use as an *in vivo* reporter of gene expression, cell lineage, protein - protein interactions and protein trafficking


The structure of GFP has been solved by several groups and was found to be an eleven-strand beta-barrel wrapped around a central alpha-helix core. This central core contains the chromophore which is spontaneously formed from a chemical reaction involving residues Ser 65, Tyr 66, and Gly 67 wherein there is cyclization of the polypeptide backbone between Ser 65 and Gly 67 to form a 5-membered ring, followed by oxidation of Tyr 66. The high quantum yield of GFP fluorescence probably arises from the nearly complete protection of the fluorophore from quenching water or oxygen molecules by burial within the beta-barrel. One of the most important attributes of GFP which makes it so useful in the life sciences is that the luminescent chromophore is formed *in vivo*, and can thus generate a labeled cellular macromolecule without the difficulties of labeling with exogenous agents.

GFP can be fused to another protein either N- or C-terminally. This is due to the fact that both termini of GFP appear rather flexible on the surface of the beta-can, so that GFP's structure is not significantly distorted or destroyed by the fused protein.



Ribbon diagram of the Green Fluorescent Protein (GFP) drawn from the wild-type crystal structure. The buried chromophore, which is responsible for GFP's luminescence, is shown in full atomic detail.

## Extrinsic fluorophores:




Normalized emission spectra of prodan in  
**1)** cyclohexane, **2)** dimethylformamide, **3)** ethanol, and  
**4)** water


The spectral sensitivity of Prodan fluorescence excitation and emission spectra to the phase state of phospholipid multilamellar vesicles has been compared in vesicles composed of pure gel, pure liquid-crystalline, and an equimolar mixture of the two phases

*Biophys J*, April 1998, p. 1984-1993, Vol. 74, No. 4


**Prodan** has both an electron-donor and an electron-acceptor substituent, resulting in a large excited-state dipole moment and extensive solvent polarity-dependent fluorescence shifts




www.probes.com



**1,8-ANS** has proved to be a sensitive probe for partially folded intermediates in protein-folding pathways. The basis of these applications is the strong fluorescence enhancement exhibited by these amphiphilic dyes when their exposure to water is lowered. Consequently, fluorescence of ANS increases substantially when proteins to which it is bound undergo transitions from unfolded to fully or partially folded states that provide shielding from water. Molten globule intermediates are characterized by particularly high ANS fluorescence intensities due to the exposure of hydrophobic core regions that are inaccessible to the dye in the native structure.




Fluorescence enhancement of 1,8-ANS (1-anilinonaphthalene-8-sulfonic acid, upon binding to protein. The image shows aqueous solutions of 1,8-ANS excited by ultraviolet light. Addition of protein (bovine serum albumin) to the solution in the cuvette on the left results in intense blue fluorescence. The fluorescence of uncomplexed free dye in the cuvette on the right is negligible in comparison.




Fluorescence emission spectra of equal concentrations of 1,8-ANS in ethanol:water mixtures. The labels adjacent to each curve indicate the percentage of ethanol in the solvent mixture.

## Nucleic acid fluorescence:


Intrinsic fluorescence is too weak and too deep in the UV for practical applications. There are numerous probes which spontaneously bind to DNA and display enhanced emission.




Ethidium Bromide  
518/605 nm




Ethidium Homodimer  
528/617 nm




Acridine Orange  
500/526 nm DNA  
460/650 nm RNA



Hoechst 33342  
350/460 nm



DAPI  
355/461 nm



TOTO  
514/533 nm

Figure 3.21. Representative DNA probes. Excitation and emission wavelengths refer to DNA-bound dye.



Ethidium bromide is one of the most widely used dyes. The mode of binding appears to be due to intercalation of the planar aromatic ring between the base pairs of double-helical DNA. EB is weakly fluorescent in water, and its intensity increases about 30-fold upon binding to DNA. Other dyes, such as acridine orange, also bind by intercalation.





Figure 9.10. Oxygen quenching of ethidium bromide in solution and when intercalated in double helical DNA. Reprinted with permission from the American Chemical Society. (From Ref. 28.)


Dye Name	Ex/Em *	Application
----------	---------	-------------

#### Dyes for Ultrasensitive Solution Quantitation

PicoGreen	502/523	Ultrasensitive reagent for solution quantitation of dsDNA
OliGreen	498/518	Ultrasensitive reagent for solution quantitation of ssDNA and oligonucleotides
RiboGreen	500/520	Ultrasensitive reagent for solution quantitation of RNA


#### Dyes for Sensitive Detection of Nucleic Acids in Gels and on Blots

SYBR Gold stain	495/537	Ultrasensitive gel stain for single- or double-stranded DNA or RNA post-electrophoresis
SYBR Green I stain	494/521	Ultrasensitive gel stain for double-stranded DNA and oligonucleotides post-electrophoresis Also useful for real-time PCR assays
SYBR Green II stain	492/513	Sensitive stain for RNA and ssDNA post-electrophoresis
SYBR DX DNA blot stain	475/499	Sensitive stain for DNA on blots (not recommended for staining RNA)



The relaxation of a single, 39  $\mu\text{m}$ -long DNA molecule stained with YOYO-1 iodide imaged at 4.5 second intervals. After the 1  $\mu\text{m}$  polystyrene sphere was trapped with optical tweezers, the attached DNA was stretched to its full extension in a fluid flow and then allowed to relax upon stoppage of fluid flow due to its entropic elasticity. The YOYO-1 iodide–DNA complex is excited with the 488 nm spectral line of the argon-ion laser and visualized through a 515 nm longpass optical filter using a Hamamatsu SIT camera with image processing.

**Perkins TT, Quake SR, Smith DE, Chu S.**  
Science 1994 May 6;264(5160):822-6



mixture of live and dead bovine sperm cells. Live sperm with intact membranes are labeled with nucleic acid stain, SYBR 14, and fluoresce green. Dead sperm, which have been killed by unprotected freeze-thawing, are labeled with propidium iodide and fluoresce red-orange.

**Molecular Probes: [www.probes.com](http://www.probes.com)**


## DNA Base Analogs:

### 2 – Aminopurine (analog of adenine)


In solution, 2-AP has a high quantum yield and a single exponential decay.


Upon incorporation into ds-DNA, its fluorescence is partially quenched and its decay becomes complex.

The sensitivity of 2-AP to its environment makes it a useful probe for studies of DNA conformation and dynamics.



### isoxanthopterin (analog of guanine)





**Principles of Fluorescence Spectroscopy**

Joseph R. Lakowicz J. Lakowicz

**Format:**Textbook Hardcover, 2nd ed., 698pp.

**ISBN:** 0306460939

**Publisher:** Kluwer Academic Publishers

**Pub. Date:** January 1999

**Topics in Fluorescence Spectroscopy** Joseph R. Lakowicz (Editor)  
(Kluwer Academic/Plenum Publishers )

**Volume 1: Techniques** January 1992

**Volume 2: Principles Volume 3: Biochemical Applications**

**Volume 4: Probe Design and Chemical Sensing**

**Volume 5: Nonlinear and Two-Photon Induced Fluorescence**

**Volume 6: Protein Fluorescence** November 2000

**Biophysical Chemistry. C. R. Cantor and P.R. Schimmel. Vol. II, (!980), Freeman, New York.**

**Applications of Fluorescence in the Biomedical Sciences. Eds. D. L. Taylor, A. S. Waggoner, F. Lanni, R. F. Murphy, and R. R. Birge. (1986), Alan R. Liss, Inc., New York.**

1. L. Stryer and R.P. Haugland. Energy Transfer: A spectroscopic ruler. Proc. Natl. Acad. Sci. USA **58**, 719 - 726 (1967).
2. L. Stryer. Fluorescence energy transfer as a spectroscopic ruler. Ann Rev. Biochem. **47**, 819 - 846 (1978).
3. Murchie, A.I. H. *et al.* Fluorescence Energy Transfer shows that the four-way DNA junction is a right -handed cross of antiparallel molecules. *Nature*, **341**, 763 - 766 (1989).
4. Tuschl, T *et al.* A three-dimensional model for the hammerhead ribozyme based on fluorescence measurements. *Science*, **266**, 785 - 789, (1994).
5. T. Heyduk and J. C. Lee. Solution studies on the structure of bent DNA in the cAMP receptor protein-lac DNA complex. *Biochemistry*, **31**, 5165 - 5171 (1992).
6. T. Heyduk and J.C. Lee. Application of fluorescence energy transfer and polarization to monitor *E. coli* cAMP receptor protein and *lac* promoter interaction. Proc. Natl. Acad. Sci. USA **87**, 1744 - 1748 (1990).
7. T. Heyduk *et al.* CAP interacts with RNA polymerase in solution in the absence of promoter DNA. *Nature*, **364**, 548 - 549 (1993)

EFTINK MR

THE USE OF FLUORESCENCE METHODS TO MONITOR UNFOLDING TRANSITIONS IN PROTEINS  
BIOPHYS J 66 (2): 482-501 FEB 1994

ROSS JBA, SCHMIDT CJ, BRAND L

TIME-RESOLVED FLUORESCENCE OF THE 2 TRYPTOPHANS IN HORSE LIVER ALCOHOL-  
DEHYDROGENASE  
BIOCHEMISTRY-US 20 (15): 4369-4377 1981